



## Toilet Dye Testing Instructions

### Things You Will Need

- Toilet cleaning products and accessories
  - Dye Testing Tablet or Colored liquid – coffee, tea, coke, juice, food coloring etc
1. Clean your toilet bowl and tank, and remove any type of sanitizing product that changes the color of your water. Flush the toilet until the water in the tank and bowl looks clear.
  2. Take the lid off the tank and set it aside in a safe place. Place dye testing tablets or colored liquid into the tank so the water achieves a deep color. Do not flush the toilet.
  3. Leave the toilet alone for 4-6 hours, and don't allow anyone to use it.
  4. Observe the water in the toilet bowl for any change in color resulting from the dyed water in the tank leaking into the bowl. If the bowl water is clear, there's no leak from the tank. If the colored liquid from the tank shows up in the bowl, most likely your flapper valve needs to be replaced.

